


	<h1>Product Specification</h1>	Specification Number	A686-15L
		Revision	01

Part No. A686-15L  
Description Lid for 15L Bucket with Spout  
Material of Construction: 316L stainless steel



Method of Construction: Handle spot welded to top of lid  
Crevice free interior


Surface Finish: Better than 1.0 microns Ra

Max. width over spout: 380mm

Diameter: 330mm

Nominal Weight: 1kg

Recommended Storage Conditions: Dry and ambient temperature

Specification Approved By:	Signature: 	Name: Kuldeep Flora Position: Quality Manager Date: 6-Aug-19
Page 1 of 1		



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

