

	<h1>Product Specification</h1>	Specification Number	A776-3000
		Revision	00

Part No. A776-3000

Description Jumbo PharmaScoop® with Ergonomic Handle

Nominal Volume: 3000ml

Material of Construction: 316L stainless steel



Method of Construction: Crevice free, all welds ground and polished

Surface Finish: Better than 0.5 microns Ra

Overall Length: 420mm


Body Length: 270mm

Body Diameter: 165mm

Nominal Weight: 1240g

BSE/TSE statement: All polishing compounds used in the manufacture of this scoop are certified BSE/TSE free

PharmaScoop is a registered trade mark belonging to Sampling Systems Ltd.

Specification Approved By:	Signature: 	Name: Kuldeep Flora Position: Quality Manager Date: 29-Jan-19
Page 1 of 1		



# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

**Wolflabs**

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

