

Sampler Specification

Specification Number

1030A

Revision

04



Part No.

1030A all lengths

Description

PowderThief - 12.5mm diameter

Part No:	Overall Length	Weight
1030A-300	370mm	350g
1030A-600	665mm	540g
1030A-1000	1065mm	795g
1030A-1500	1565mm	1020g
1030A-2000	2065mm	1485g



Materials of Construction:

Inner Rod:	316L stainless steel
Handles:	316L stainless steel
Handle Boss:	316L stainless steel
Outer Tube:	316L stainless steel
Clutch Housing:	316L stainless steel
Clutch Socket:	316L stainless steel
Clutch:	PTFE

External Polish:

<1 microns Ra

Internal Polish:


None

Recommended Storage Conditions:

Dry and ambient temperature

Spare Parts:

PTFE Clutch	Part No: SS002300
Clutch Housing	Part No: SS001900

Specification Approved By:	Signature: 	Name: Kuldeep Flora Position: Quality Manager Date: 17-Feb-20
Page 1 of 1		



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.