

Sampler Specification

Specification Number

1030B



Revision

04

Part No. 1030B all lengths

Description PowderThief - 19mm diameter

Part No:	Overall Length	Weight
1030B-600	680mm	1075g
1030B-1000	1080mm	1580g
1030B-1500	1580mm	2080g
1030B-2000	2080mm	2620g



Materials of Construction:

Inner Rod: 316L stainless steel
 Handles: 316L stainless steel
 Handle Boss: 316L stainless steel
 Outer Tube: 316L stainless steel
 Clutch Housing: 316L stainless steel
 Clutch Socket: 316L stainless steel
 Clutch: PTFE
 Disk: PTFE

External Polish: <1 microns Ra

Internal Polish: None

Recommended Storage Conditions: Dry and ambient temperature

Spare Parts: PTFE Clutch Part No: SS002301
 Clutch Housing Part No: SS001901
 PTFE Disk Part No: SS002001

Specification Approved By:	Signature: 	Name: Kuldeep Flora Position: Quality Manager Date: 17-Feb-20
Page 1 of 1		



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.