

Sampler Specification

Specification Number

1030C



Revision

04

Part No. 1030C all lengths

Description PowderThief - 25mm diameter

Part No:	Overall Length	Weight
1030C-600	680mm	1250g
1030C-1000	1080mm	1780g
1030C-1500	1580mm	2510g


Materials of Construction:


Inner Rod:	316 stainless steel
Handles:	316 stainless steel
Handle Boss:	316 stainless steel
Outer Tube:	316 stainless steel
Clutch Housing:	316 stainless steel
Clutch Socket:	316 stainless steel
Clutch:	PTFE
Disk:	PTFE

External Polish: <1 microns Ra

Internal Polish: None

Recommended Storage Conditions: Dry and ambient temperature

Spare Parts:	PTFE Clutch	Part No: SS002301
	Clutch Housing	Part No: SS001902
	PTFE Disk	Part No: SS002002

Specification Approved By:	Signature: 	Name: Kuldeep Flora Position: Quality Manager Date: 17-Feb-20
Page 1 of 1		



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.