


1. Part Number: 1210A-1000

2. Photo:



3. Product Information:

Description	Slot Sampler 1000mm
Materials:	316L stainless steel
External Polish:	Better than 0.8 microns Ra
Internal Polish:	None
Outer Diameter of Outer Tube:	25mm
Inner Diameter of Inner Tube:	17mm
Overall Length:	1025mm
Slot Length:	120mm
Slot Width:	15mm
Number of Slots	4
Slot Positions:	



Nominal Weight of Sampler:	1.9kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 20-Mar-20



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

