


**1. Part Number:** 1220J-1500

**2. Photo:**



**3. Product Information:**

Description	Jumbo Pocket Sampler 1500mm
Materials:	316L stainless steel
External Polish:	Better than 0.8 microns Ra
Internal Polish:	None
Overall Length:	1570mm
Diameter of Outer Tube:	40mm
Pocket Length:	125mm
Pocket Width:	20mm
Nominal Pocket Volume:	90ml
Number of Pockets:	4
Slot Positions:	



Nominal Weight:	5.6 kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 25-Mar-20



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

