

1. Part Number: 1320A-01

2. Photo:



3. Product Information:

Description	Sack Master – with Plain End
Method of Construction:	Single piece construction
Materials:	316L stainless steel
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	400mm
Outer Tube Diameter:	25mm
Inner Diameter:	22mm
Nominal Weight:	335g
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 4-Oct-19



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.