

1. Part Number: 1640A-1500

2. Photograph:



3. Product Information:

Description:	Point Sampler 1500mm
Material of Construction:	316L stainless steel
Method of Construction:	Welds ground and polished
Finish	- External Bright Polish - Internal Machine finish
Nominal Sample Volume:	100 ml
Overall Length:	1730 mm
Max. Diameter:	42 mm (NB. The scraper will add an extra 10mm)
Nominal Weight of Sampler:	2.1 Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this item are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Nov-20



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

