

- 1. Part Number:** 2000A-1000
2000T-025
2000T-050
2000T-075
2000T-100

2. Photo:



3. Product Information:

Description: U-D Cohesive Sampler & Tips
Materials: 316L stainless steel
(Note: the stainless steel spring is not grade 316L)
Method of Construction: Crevice free, all welds ground and polished

Sampler Body (Part No: 2000A-1000)

Max. Diameter: 26mm
Overall Length (sampler & tip): 1080mm
External Finish: <1 Micron Ra
Nominal Weight of Sampler (2000A-1000): 660g

Tips

Part No.	Length	Overall Diameter	Depth of Sample Chamber	Weight
2000T-025	30mm	26mm	3.2mm	60g
2000T-050	30mm	26mm	5.5mm	57g
2000T-075	30mm	26mm	7.6mm	56g
2000T-100	35mm	26mm	9.8mm	74g

BSE/TSE statement: All polishing compounds used in the manufacture of this sampler are of vegetable origin and so are BSE/TSE free.

Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 7-Oct-20



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

