


**1. Part Number:** 2015A-10

**2. Photo:**



### 3. Product Information:

|                                 |  |
|---------------------------------|--|
| Description                     | Cohesive Sampler 10ml  |
| Materials:                      | 316L stainless steel   |
| Method of Construction:         | Crevice free, all welds ground and polished  |
| Surface Finish:                 | Better than 0.8 microns Ra   |
| Nominal Sample Volume:          | 10ml   |
| Overall Length:                 | 1030mm   |
| Diameter of Outer Tube:         | 25mm   |
| Scraper Length                  | 110mm  |
| Nominal Weight of Sampler:      | 1.9Kg  |
| BSE/TSE statement:              | All polishing compounds used in the manufacture of this sampler are of vegetable origin and so are BSE/TSE free. |
| Recommended Storage Conditions: | Dry and ambient temperature  |

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Oct-20



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

