


**1. Part Number:** 2210A-1000

**2. Photo:**



### 3. Product Information:

Description	Cohesive Slot Sampler 1000mm
Materials:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Outer Diameter of Outer Tube:	25mm (plus 10mm approx. for scrapers)
Inner Diameter of Inner Tube:	17mm
Overall Length:	1025mm
Slot Length:	120mm
Slot Width:	15mm
Number of Slots	4
Slot Positions:	



External Polish:	< 0.8 Microns Ra
Internal Polish:	None
Nominal Weight of Sampler:	2.0Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin and so are BSE/TSE free.
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Oct-20



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.