

1. Part Number: 5035T-1501

2. Photo:



3. Product Information:

Description	Liquid Master PTFE – 1500mm
Materials of Construction:	PTFE & FEP (for product contact parts) <ul style="list-style-type: none"> ○ FDA 21 CFR 177.1550 ○ USP Class VI ○ EU regulation 10/2011 ○ EC regulation 1935/2004
Outer Diameter of Outer Tube:	25mm
Internal Diameter of Outer Tube:	19mm
Max. Diameter of Tip:	32mm
Diameter of hole in end of tip:	10mm
Overall Length:	1500mm
BSE/TSE Free:	Yes
Nominal Weight:	860g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 4-Oct-19



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

