

	<h1>Product Specification</h1>	Specification Number	5055A-50S
		Revision	03

Part No. 5055A-50S

Description Disposable Thumb Tubes

Material of Construction: High Density Polyethylene (HDPE)

- FDA approved CFR 177.1520
- Conforms to EU regs 10/2011



Method of Construction: Single piece, Extrusion

Overall Diameter: 12mm

Length: 1000mm


Nominal Weight (each): 44g

Individually Bagged? Yes (heat sealed PE bag)

Method of Sterilisation: Gamma Irradiation (25-45kGy)

Material BSE/TSE Free Yes

Number of Samplers per Box: 50

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 20-Feb-13
Page 1 of 1		



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.