

**1. Part Number:** 8284B-150S

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® Cupped Spatula Blue – 150mm - Sterile
Use:	Single use, disposable spatula
Material of Construction:	Polypropylene (PP) – rigid, blue <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulation 10/2011</li> <li>• Conforms to EC Regulation 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Blade Shape:	U
Blade Style:	Closed Tip
Nominal Volume (flat filled):	21ml
Nominal Length:	240mm
Blade Length:	155mm
Blade Width:	18mm
Nominal Weight:	18g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25 to 45 kGy)
Quantity per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This item has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable device. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 26-Oct-20





**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

