

Sartorius MA Infrared Moisture Analysers Continued

SARTORIUS

Accessories

Code	Alt Ref	Description	Price
BAL5040	YSS43-02	External calibration weight	£190.00
BAL5042	YTM15MA	Temperature adjustment set	£961.00
BAL4730	YDP30	Thermo transfer and thermo direct printer	£781.00
BAL4732	69Y03287	Standard thermal paper for YDP30	£33.40

OHAUS MB90 and MB120 Moisture Analysers



Geared to increase precision and efficiency in the lab, these Moisture Analysers feature a precisely controlled halogen heating system which helps to yield fast and repeatable results with readability up to 0.01%. The robust housing and die-cast structure is designed for easy, tool-free cleaning and durability.

The MB90 is ideal for routine moisture analysis applications. It begins to dry samples in seconds and performs up to 40% faster than traditional infrared methods, providing fast, precise results that you can rely on. So easy to operate and clean - you can test more samples in less time!

The MB120 is designed for professional moisture analysis with advanced features including a temperature guide which analyses a sample to determine optimum drying temperature and a Methods Library that can store up to 100 methods for faster testing. This high-performance moisture analyser providing precise, repeatable results of 0.01% readability can streamline your daily workflow.



Code	Alt Ref	Model	Capacity, g	Temp Range, °C	Readability, mg	Pan Dims, mm	Moisture Range, %	Price
BAL5532	30303303	MB90	0.5 - 90	40-200	1/0.01%	90	3 - 99	£2115.00
BAL5534	30303304	MB120	0.5 - 120	40-230	1/0.01%	90	4 - 99	£3085.00

Accessories

Code	Alt Ref	Description	Pack	Price
BAL5536	30246947	Pan holder	Each	£36.00
BAL5538	30284478	In-use cover	Each	£31.00
BAL5540	30284477	Scoop	Each	£25.00
BAL6240	11113857	Temperature calibration kit	Each	£282.00
BAL5152	80850043	Security device	Each	£138.00
BAL4234	80850086	Aluminum pans	80	£29.00
BAL6236	80850087	Glass fiber filters	200	£44.00
BAL5542	80252477	Sample cage	Each	£36.00
BAL5544	80252478	Reusable pans, 7mm	3	£76.00
BAL5546	80252479	Reusable pans, 14mm	3	£78.00
BAL1021	30064202	Printer SF40A	Each	£418.00



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.