

SciSpin

Centrifuges

SciSpin MINI Microfuges



SciSpin MINI

- Ultra compact design
- Four lid colour options: green, blue, yellow & pink
- Max 7,000 RPM / 2,680 g
- Includes 8 x 1.5 / 2.2 ml & 16 x PCR tubes
- Dual lid lock for safety
- Whisper quiet, <45dB
- Ideal for quick spin downs and micro filtration
- Rapid acceleration and deceleration
- Close lid to begin run & open to stop

Ordering Information

SciSpin MINI Microfuge

Part No.	Description
SQ-6050	SciSpin MINI Microfuge blue lid, including two rotors of SQ-6058 & SQ-6059
SQ-6051	SciSpin MINI Microfuge green lid, including two rotors of SQ-6058 & SQ-6059
SQ-6052	SciSpin MINI Microfuge pink lid, including two rotors of SQ-6058 & SQ-6059
SQ-6053	SciSpin MINI Microfuge yellow lid, including two rotors of SQ-6058 & SQ-6059

Accessories

Part No.	Description
SQ-6058	Rotor for 8 x 1.5 / 2.2 ml, (7,000 RPM / 2,680 g)
SQ-6059	Rotor for 16 0.2 ml PCR tubes, 2 x 8 strips (7,000 RPM / 2,680 g)

SciSpin MINI

Rotor Capacity	8 x 1.5 / 2.2 ml
	16 x 0.2ml PCR
Max Speed	7,000 RPM
Max RCF	2,680 g
Timer	Continuous operation
Dimensions (WxDxH)	150mm x 150mm x 117mm
Weight	0.5 kg
Noise level	45 dB
Power	100 - 240V 50/60Hz
Noise level	56dB





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.