

# SciSpin

Centrifuges

## SciSpin ONE Centrifuge



## SciSpin ONE

- Specifically designed for small batch routine applications
- Capacity: 12 x vacutainers/blood tubes  
12 x 10ml tubes  
8 x 15ml conical tubes
- Compact design
- Speed Range: 300 – 4,500 RPM
- Timer from 30 sec - 99 min or continuous run
- Pulse button for quick spins
- Brushless motor
- Rapid acceleration and gentle braking
- Whisper quiet, <56dB
- Over speed and over temperature detection

### Ordering Information

#### SciSpin ONE Clinical Centrifuge

Part No.	Description
SS-6100	SciSpin ONE Clinical Centrifuge with rotor SS-6105

#### Accessories

Part No.	Description
SS-6105	Plastic rotor, 4500rpm, 15ml culture tubes x8
SS-6108	Rotor adapter, for 4 - 7 ml tubes used with SS-6105, 12pcs/pk
SS-6109	Rotor adapter Mat, used with SS-6105, 12pcs/pk

#### SciSpin ONE

Rotor Capacity	8 x 15ml conical tubes 12 x 10ml, 7ml or 5ml tubes (inc. blood tubes)
Max Speed	4,500 RPM
Max RCF	2,490 g
Timer	30 sec - 99 min / continuous run
Dimensions (WxDxH)	255 mm x 245 mm x 140 mm
Weight	6 kg
Noise level	<56 dB
Power	100 - 240V 50/60Hz



SS-6108 Adapter



SS-6109 Adapter mat



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.