




## Azlon® Hinged safety shield

### For **excellent protection** in hazardous situations

Incredibly strong and lightweight, the Azlon® hinged type safety shield can be clamped to the bench for increased stability.

- Excellent protection in potentially hazardous situations
- Adjustable - can be used as a tight triangle or extended out to a straight line configuration
- Incredibly strong and lightweight
- Supplied with two pairs of plastic coated metal feet
- Feet may be clamped or bolted onto the bench to give greater stability



### Product Specification

Item Code	Centre panel width mm	Height cm	Outside panel widths mm	Total width mm	Pack quantity
SGT002	230	38	180	590	1
SGT004	230	61	180	590	1
SGT006	300	61	226	750	1

