

Sterilin Microtiter Plates, Black and White

Sterilin black 96-well plates are for fluorescence work, and white plates are for bioluminescence studies



Product highlights

- SBS

Details

- Black 96-well plates are ideal for fluorescence work where the black pigment prevents 'cross-talk' between the wells — the low background fluorescence minimizes light scattering
- White 96-well plates are ideal for bioluminescence studies where the pigment aids reflectivity and sensitivity
- SBS compliant, ensuring suitable fit with all automation
- Manufactured from virgin polystyrene
- Supplied non-sterile, providing a medium binding surface
- Chimney well design to help reduce contamination



Note: Supplied without a lid; if required, the 642000 lid is a suitable fit.

Sterilin Microtiter Plates, Black and White

| Description | Well capacity | Material | Lid | Sterility | Inner pack qty | Case Qty. | Cat. No. |
|----------------------------------|---------------|----------|-----|-----------|----------------|-----------|----------|
| Black 96-well plate, Flat bottom | 400µL | PS | No | NS | 5 | 50 | 611F96BK |
| White 96-well plate, Flat bottom | 400µL | PS | No | NS | 5 | 50 | 611F96WT |



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.