


Thermo Scientific™ Nunc™ Immuno MicroWell™ Plates

Thermo
SCIENTIFIC

- Plates intended for quantitative and qualitative solid phase immuno assays (e.g. ELISA) and binding assays
- Choice of flat, round bottom and flat with curved edges for optimal washing and liquid retrieval
- Raised rims on the wells
- Certified reproducibility of binding
- Alphanumeric orientation system
- Without lid, non-sterile

Code	Surface Treatment	Colour	Description	Well Volume, μL	Bottom	Inner Pack	Pack	Price
475094	PolySorp	Clear	Pinchbar	400	Flat	5	60	£154.00
456529	PolySorp	Clear	High flange	400	Flat	10	180	£422.00
446140	PolySorp	Clear	Pinchbar	350	C	5	60	£144.00
475434	PolySorp	Clear	High flange	300	U	5	60	£147.00
467320	MediSorp	Clear	High flange	400	Flat	5	60	£213.00
467340	MultSorp	Clear	High flange	400	Flat	5	60	£213.00
456537	MaxiSorp	Clear	High flange	400	Flat	10	180	£566.00
439454	MaxiSorp	Clear	Pinchbar	400	Flat	5	60	£210.00
442404	MaxiSorp	Clear	Pinchbar	400	Flat	5	60	£155.00
430341	MaxiSorp	Clear	Pinchbar	350	C	5	60	£144.00
446612	MaxiSorp	Clear	Pinchbar	350	C	5	60	£212.00
449824	MaxiSorp	Clear	High flange	300	U	5	60	£142.00
436110	MaxiSorp	White	Pinchbar	400	Flat	10	80	£334.00
437111	MaxiSorp	Black	Pinchbar	400	Flat	10	80	£332.00



Thermo Scientific™ Plates with Covalent or Affinity Binding Surfaces

Thermo
SCIENTIFIC

- Increase the sensitivity of your ELISA techniques by using either covalent or affinity capture binding to immobilise your molecule to the microtitre plate
- Covalent Binding Surfaces chemically attach the molecule to the plate from the NH₂ end of the molecule, this means all molecules are oriented the same way and are impossible to wash away
- Affinity Binding Surfaces using streptavidin allow biotinylated molecules to bind strongly to the surface also in the same orientation
- Both methods do not require a blocking step for the ELISA
- No blocking required speeding up assay time
- Significantly reduces the amount of molecule required for coating
- Well volume: 400 μL
- 5 year shelf life of a pre-coated plate

Code	Surface Treatment	Colour	Pack	Price
436006	Amine	Clear	30	£613.00
436007	Amine	White	30	£885.00
436008	Amine	Black	30	£885.00
436015	Streptavidin	White	15	£594.00
436016	Streptavidin	Black	15	£594.00
436014	Streptavidin	Clear	15	£456.00



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.