

thermo scientific

Thermo Scientific ClipTip Pipetting System

Sealed in security. Feel the difference.

Selection Guide



ThermoFisher
SCIENTIFIC

Worry-free pipetting is finally here

Eliminate the frustration of banging tips on your pipette


With a light touch, each tip is clipped securely on the pipette.

Reduce the risk of loose or leaking tips

When Thermo Scientific™ ClipTip™ Pipette Tips are on, they are sealed and will not loosen. The ClipTip Pipetting System remains sealed until you decide to eject tips, even in the most comprehensive applications, including multiple tip touches and mixing cycles. This enables consistent and reproducible pipetting results from user to user.

ClipTip Pipette Tips

- Tips that clip securely on the pipette.
- Versatile product range with 12.5-1250 µL 96-format tips and 12.5-125µl 384-format tips.
- Several packaging options from racked to environmentally friendly and space saving reload towers.
- All ClipTip tips are certified free of RNase, DNase, DNA, ATP and endotoxin contamination.
- Low retention technology minimizes liquid retention, and enable maximum sample recovery. Available in ClipTip 20-1250 96-format tips.
- Convenient color coding makes it easy to find the compatible F1-ClipTip or E1-ClipTip pipettes.
- Extended length tips available in CT12.5 Ext and CT300 Ext sizes to help enable better reach with specific labware formats.



F1-ClipTip™ Manual Pipettes

- Interlocking ClipTip technology provides a secure seal on every channel.
- Extremely light tip attachment and ejection forces to help reduce strain injuries.
- Effortless and secure volume adjustment with light and soundless volume adjustment mechanism and volume lock.
- Light-weight handle design with adjustable finger rest.
- Complete range of single, 8- and 12-channel pipettes within volume range of 0.1-1000 µL.

E1-ClipTip™ Electronic Pipettes

- Interlocking ClipTip technology provides a secure seal on every channel.
- Adjustable tip-spacing combined with the electronic tip ejection, and index finger operation ensures more ergonomic and efficient pipetting experience.
- An intuitive personalized user-interface makes even the most complex protocols easier, and efficient to complete.
- Versatile product range from single, 8-, 12- and 16-channel to 6-12-channel adjustable tip spacing pipettes.

F1-ClipTip manual pipettes and ClipTip pipette tips ordering information

| Color key/description | Volume range, μL | Increment/Packaging | Cat. No. |
|---|-----------------------------|---|----------|
| 2/10 μL Pipettes | | | |
| Increment | | | |
| • F1-ClipTip 0.1-2 μL | 0.1-2 μL | 0.002 μL | 4641310N |
| • F1-ClipTip 1-10 μL | 1-10 μL | 0.02 μL | 4641320N |
| • F1-ClipTip 8-ch 1-10 μL | 1-10 μL | 0.02 μL | 4661210N |
| • F1-ClipTip 12-ch 1-10 μL | 1-10 μL | 0.02 μL | 4661220N |
| 12.5 Ext μL Tips* | | | |
| Packaging | | | |
| • ClipTip 12.5 Ext | 0.1-12.5 μL | Racked, 96 tips/rack, 10 racks/pack | 94410060 |
| • ClipTip 12.5 Ext | 0.1-12.5 μL | Racked, Sterile, 96 tips/rack, 10 racks/pack | 94410063 |
| • ClipTip 12.5 Ext | 0.1-12.5 μL | Racked, Filtered, Sterile, 96 tips/rack, 10 racks/pack | 94420063 |
| 20 μL Pipettes | | | |
| Increment | | | |
| • F1-ClipTip 2-20 μL | 2-20 μL | 0.02 μL | 4641180N |
| 20 μL Tips[†] | | | |
| Packaging | | | |
| • ClipTip 20 | 1-20 | Racked, 96 tips/rack, 10 racks/pack | 94410210 |
| • ClipTip 20 | 1-20 | Racked, sterile, 96 tips/rack, 10 racks/pack | 94410213 |
| • ClipTip 20 | 1-20 | Reload tower, 96 tips/insert, 10 inserts/pack | 94410217 |
| • ClipTip 20 | 1-20 | Reload tower, sterile, 96 tips/insert, 10 inserts/pack | 94410218 |
| • ClipTip 20 | 1-20 | Racked, filtered, sterile, 96 tips/rack, 10 racks/pack | 94420213 |
| • ClipTip 20 | 1-20 | Reload insert, filtered, sterile, 96 tips/insert, 10 inserts/pack | 94420218 |
| 50 μL Pipettes | | | |
| Increment | | | |
| • F1-ClipTip 5-50 μL | 5-50 | 0.1 μL | 4641190N |
| • F1-ClipTip 8-ch 5-50 μL | 5-50 | 0.1 μL | 4661120N |
| • F1-ClipTip 12-ch 5-50 μL | 5-50 | 0.1 μL | 4661160N |
| 50 μL Tips | | | |
| Packaging | | | |
| • ClipTip 50 | 5-50 | Racked, 96 tips/rack, 10 racks/pack | 94410250 |
| • ClipTip 50 | 5-50 | Racked, sterile, 96 tips/rack, 10 racks/pack | 94410253 |
| • ClipTip 50 | 5-50 | Reload tower, 96 tips/insert, 10 inserts/pack | 94410257 |
| • ClipTip 50 | 5-50 | Reload tower, sterile, 96 tips/insert, 10 inserts/pack | 94410258 |
| • ClipTip 50 | 5-50 | Racked, filtered, sterile, 96 tips/rack, 10 racks/pack | 94420253 |
| • ClipTip 50 | 5-50 | Reload insert, filtered, sterile, 96 tips/insert, 10 inserts/pack | 94420258 |
| 100 μL / 200 μL Pipettes | | | |
| Increment | | | |
| • F1-ClipTip 10-100 μL | 10-100 | 0.2 μL | 4641200N |
| • F1-ClipTip 20-200 μL | 20-200 | 0.2 μL | 4641210N |
| • F1-ClipTip 8-ch 10-100 μL | 10-100 | 0.2 μL | 4661130N |
| • F1-ClipTip 12-ch 10-100 μL | 10-100 | 0.2 μL | 4661170N |
| 200 μL Tips | | | |
| Packaging | | | |
| • ClipTip 200 | 10-200 | Racked, 96 tips/rack, 10 racks/pack | 94410310 |
| • ClipTip 200 | 10-200 | Racked, sterile, 96 tips/rack, 10 racks/pack | 94410313 |
| • ClipTip 200 | 10-200 | Reload tower, 96 tips/insert, 10 inserts/pack | 94410317 |
| • ClipTip 200 | 10-200 | Reload tower, sterile, 96 tips/insert, 10 inserts/pack | 94410318 |
| • ClipTip 200 | 10-200 | Racked, filtered, sterile, 96 tips/rack, 10 racks/pack | 94420313 |
| • ClipTip 200 | 10-200 | Reload insert, filtered, sterile, 96 tips/insert, 10 inserts/pack | 94420318 |

* CT 12.5 Ext tips and compatible F1-ClipTip pipettes utilize 'snap and seal' technology. All ClipTip pipette tips are low retention except for the CT12.5, CT12.5 Ext and 384-format tips.


[†] CT20 tips are also used with old F1-ClipTip 1-10 μL pipette models with catalogue numbers 4641170N, 4661110N, 4661150N and 4651190N.

| | | |
|--|---|--|
| ● 10 μL / 12.5 μL / 20 μL | ● 100 μL / 125 μL / 200 μL | ● 850 μL / 1000 μL |
| ● 30 μL / 50 μL | ● 300 μL | ● 1250 μL |

F1-ClipTip Manual Pipettes and ClipTip Pipette Tips ordering information

| Color key/description | Volume range, μL | Increment/Packaging | Cat. No. |
|--|-----------------------------|---|----------|
| 300 μL Pipettes | | | |
| | | Increment | |
| • F1-ClipTip 30-300 μL | 30-300 | 1 μL | 4641220N |
| • F1-ClipTip 8-ch 30-300 μL | 30-300 | 1 μL | 4661140N |
| • F1-ClipTip 12-ch 30-300 μL | 30-300 | 1 μL | 4661180N |
| | | Packaging | |
| • ClipTip 300 | 30-300 | Racked, 96 tips/rack, 10 racks/pack | 94410510 |
| • ClipTip 300 | 30-300 | Racked, sterile, 96 tips/rack, 10 racks/pack | 94410513 |
| • ClipTip 300 | 30-300 | Reload tower, 96 tips/insert, 10 inserts/pack | 94410517 |
| • ClipTip 300 | 30-300 | Reload tower, sterile, 96 tips/insert, 10 inserts/pack | 94410518 |
| • ClipTip 300 | 30-300 | Racked, filtered, sterile, 96 tips/rack, 10 racks/pack | 94420513 |
| • ClipTip 300 | 30-300 | Reload insert, filtered, sterile, 96 tips/insert, 10 inserts/pack | 94420518 |
| • ClipTip 300 Ext | 30-300 | Racked, 96 tips/rack, 8 racks/pack | 94410610 |
| • ClipTip 300 Ext | 30-300 | Racked, sterile, 96 tips/rack, 8 racks/pack | 94410613 |
| • ClipTip 300 Ext | 30-300 | Reload tower, 96 tips/insert, 10 inserts/pack | 94410617 |
| • ClipTip 300 Ext | 30-300 | Racked, filtered, sterile, 96 tips/rack, 8 racks/pack | 94420613 |
| | | Increment | |
| • F1-ClipTip 100-1000 μL ** | 100-1000 | 1 μL | 4641230N |
| | | Packaging | |
| • ClipTip 1000 | 100-1000 | Racked, 96 tips/rack, 8 racks/pack | 94410710 |
| • ClipTip 1000 | 100-1000 | Racked, sterile, 96 tips/rack, 8 racks/pack | 94410713 |
| • ClipTip 1000 | 100-1000 | Reload tower, 96 tips/insert, 8 inserts/pack | 94410717 |
| • ClipTip 1000 | 100-1000 | Reload tower, sterile, 96 tips/insert, 8 inserts/pack | 94410718 |
| • ClipTip 1000 | 100-1000 | Racked, filtered, sterile, 96 tips/rack, 8 racks/pack | 94420713 |
| • ClipTip 1000 | 100-1000 | Reload insert, filtered, sterile, 96 tips/insert, 8 inserts/pack | 94420718 |

**Short Shaft



F1-ClipTip manual fixed volume pipettes ordering information

| Description | Volume Range, μL | Cat. No. |
|----------------------------------|-----------------------------|----------|
| F1-ClipTip 10 μL | 10 μL | 4651280N |
| F1-ClipTip 20 μL | 20 μL | 4651200N |
| F1-ClipTip 25 μL | 25 μL | 4651210N |
| F1-ClipTip 50 μL | 50 μL | 4651220N |
| F1-ClipTip 100 μL | 100 μL | 4651230N |
| F1-ClipTip 200 μL | 200 μL | 4651240N |
| F1-ClipTip 250 μL | 250 μL | 4651250N |
| F1-ClipTip 500 μL ** | 500 μL | 4651260N |
| F1-ClipTip 1000 μL ** | 1000 μL | 4651270N |

**Short Shaft



F1-ClipTip GLP Kit 2

F1-ClipTip Good Laboratory Pipetting (GLP) kits ordering information

| Description | Contents of Kit | Cat. No. |
|-----------------------------|---|----------|
| F1-ClipTip GLP Kit 1 | <ul style="list-style-type: none"> • F1-ClipTip 1-10 μL single channel • F1-ClipTip 10-100 μL single channel • F1-ClipTip 100-1000 μL single channel • F-Series Stand • ClipTip 12.5 Ext rack • ClipTip 200 rack • ClipTip 1000 rack • 25 mL reservoir • 25 mL divided reservoir • 100 mL divided reservoir | 4701140N |
| F1-ClipTip GLP Kit 2 | <ul style="list-style-type: none"> • F1-ClipTip 2-20 μL single channel • F1-ClipTip 20-200 μL single channel • F1-ClipTip 100-1000 μL single channel • F1-ClipTip 30-300 μL 8-channel • F-Series Stand • ClipTip 20 rack • ClipTip 200 rack • ClipTip 300 rack • ClipTip 1000 rack • 25 mL reservoir • 25 mL divided reservoir • 100 mL divided reservoir | 4701150N |

E1-ClipTip Electronic Pipettes** and ClipTip Pipette Tips ordering information

| Color Key/Description | Volume Range, μL | Tip Spacing Range, mm | Increment/Packaging | Cat. No. |
|---|-----------------------------|-----------------------|---|-----------|
| 12.5 μL Pipettes* | | | | |
| Increment | | | | |
| • E1-ClipTip 0.5-12.5 μL | 0.5-12.5 | | 0.01 μL | 467000BT |
| • E1-ClipTip 8-ch 0.5-12.5 μL | 0.5-12.5 | | 0.01 μL | 4671000BT |
| • E1-ClipTip 12-ch 0.5-12.5 μL | 0.5-12.5 | | 0.01 μL | 4671010BT |
| 12.5 μL Tips* | | | | |
| Packaging | | | | |
| • ClipTip 12.5 | 0.5-12.5 | | Racked, 96 tips/rack, 10 racks/pack | 94410040 |
| • ClipTip 12.5 | 0.5-12.5 | | Racked, sterile, 96 tips/rack, 10 racks/pack | 94410043 |
| • ClipTip 12.5 | 0.5-12.5 | | Racked, filtered, sterile, 96 tips/rack, 10 racks/pack | 94420043 |
| • ClipTip 12.5 Ext | 0.5-12.5 | | Racked, 96 tips/rack, 10 racks/pack | 94410060 |
| • ClipTip 12.5 Ext | 0.5-12.5 | | Racked, sterile, 96 tips/rack, 10 racks/pack | 94410063 |
| • ClipTip 12.5 Ext | 0.5-12.5 | | Racked, filtered, sterile, 96 tips/rack, 10 racks/pack | 94420063 |
| 384, 12.5 μL Pipettes* | | | | |
| Increment | | | | |
| • E1-ClipTip 16-ch 0.5-12.5 μL | 0.5-12.5 | | 0.01 μL | 4671020BT |
| • E1-ClipTip Equalizer 384 8-ch 0.5-12.5 μL | 0.5-12.5 | 4.5-14.2 | 0.01 μL | 4672010BT |
| • E1-ClipTip Equalizer 384 12-ch 0.5-12.5 μL | 0.5-12.5 | 4.5-9.0 | 0.01 μL | 4672020BT |
| 384, 12.5 Tips* | | | | |
| Packaging | | | | |
| • ClipTip 384 12.5 | 0.5-12.5 | | Racked, 384 tips/rack, 10 racks/pack | 94410050 |
| • ClipTip 384 12.5 | 0.5-12.5 | | Racked, sterile, 384 tips/rack, 10 racks/pack | 94410053 |
| • ClipTip 384 12.5 | 0.5-12.5 | | Racked, filtered, sterile, 384 tips/rack, 10 racks/pack | 94420053 |
| • ClipTip 384 12.5 Ext | 0.5-12.5 | | Racked, 384 tips/rack, 10 racks/pack | 94410070 |
| • ClipTip 384 12.5 Ext | 0.5-12.5 | | Racked, sterile, 384 tips/rack, 10 racks/pack | 94410073 |
| • ClipTip 384 12.5 Ext | 0.5-12.5 | | Racked, filtered, sterile, 384 tips/rack, 10 racks/pack | 94420073 |
| 384, 30 μL Pipettes* | | | | |
| Increment | | | | |
| • E1-ClipTip 16-ch 1-30 μL | 1-30 | | 0.01 μL | 4671030BT |
| • E1-ClipTip Equalizer 384 8-ch, 1-30 μL | 1-30 | 4.5-14.2 | 0.01 μL | 4672030BT |
| • E1-ClipTip Equalizer 384 12-ch, 1-30 μL | 1-30 | 4.5-9.0 | 0.01 μL | 4672040BT |
| 384, 30 μL Tips* | | | | |
| Packaging | | | | |
| • ClipTip 384 30 | 1-30 | | Racked, 384 tips/rack, 10 racks/pack | 94410100 |
| • ClipTip 384 30 | 1-30 | | Racked, sterile, 384 tips/rack, 10 racks/pack | 94410103 |
| • ClipTip 384 30 | 1-30 | | Racked, filtered, sterile, 384 tips/rack, 10 racks/pack | 94420103 |
| 125 μL Pipettes | | | | |
| Increment | | | | |
| • E1-ClipTip 1-ch 2-125 μL | 2-125 | | 0.1 μL | 4670020BT |
| • E1-ClipTip 8-ch 2-125 μL | 2-125 | | 0.1 μL | 4671040BT |
| • E1-ClipTip 12-ch 2-125 μL | 2-125 | | 0.1 μL | 4671050BT |
| • E1-ClipTip Equalizer 8-ch 2-125 μL | 2-125 | 9.0-14.2 | 0.1 μL | 4672050BT |
| 200 μL Tips | | | | |
| Packaging | | | | |
| • ClipTip 200 | 2-200 | | Racked, 96 tips/rack, 10 racks/pack | 94410310 |
| • ClipTip 200 | 2-200 | | Racked, sterile, 96 tips/rack, 10 racks/pack | 94410313 |
| • ClipTip 200 | 2-200 | | Reload tower, 96 tips/insert, 10 inserts/pack | 94410317 |
| • ClipTip 200 | 2-200 | | Reload tower, sterile, 96 tips/insert, 10 inserts/pack | 94410318 |
| • ClipTip 200 | 2-200 | | Racked, filtered, sterile, 96 tips/rack, 10 racks/pack | 94420313 |
| • ClipTip 200 | 2-200 | | Reload insert, filtered, sterile, 96 tips/insert, 10 inserts/pack | 94420318 |
| 384, 125 μL Pipettes* | | | | |
| Increment | | | | |
| • E1-ClipTip 16-ch 2-125 μL | 2-125 | | 0.1 μL | 4671060BT |
| • E1-ClipTip Equalizer 384 8-ch 2-125 μL | 2-125 | 4.5-14.2 | 0.1 μL | 4672060BT |
| • E1-ClipTip Equalizer 384 12-ch 2-125 μL | 2-125 | 4.5-9.0 | 0.1 μL | 4672070BT |
| 384, 125 μL TIPS* | | | | |
| Packaging | | | | |
| • ClipTip 384 125 | 2-125 | | Racked, 384 tips/rack, 10 racks/pack | 94410150 |
| • ClipTip 384 125 | 2-125 | | Racked, sterile, 384 tips/rack, 10 racks/pack | 94410153 |
| • ClipTip 384 125 | 2-125 | | Racked, filtered, sterile, 384 tips/rack, 10 racks/pack | 94420153 |

● 2 μL /10 μL / 12.5 μL / 20 μL
● 100 μL / 125 μL / 200 μL
● 850 μL / 1000 μL
● 30 μL / 50 μL
● 300 μL
● 1250 μL

E1-ClipTip Electronic Pipettes** and ClipTip Pipette Tips ordering information

| Color Key/Description | Volume Range, μL | Tip Spacing Range, mm | Increment/Packaging | Cat. No. |
|---|-----------------------------|-----------------------|---|-----------|
| 300 μL Pipettes | | | | |
| | | | Increment | |
| • E1-ClipTip 1-ch 10-300 μL | 10-300 | | 0.1 μL | 4670030BT |
| • E1-ClipTip 8-ch 10-300 μL | 10-300 | | 0.1 μL | 4671070BT |
| • E1-ClipTip 12-ch 10-300 μL | 10-300 | | 0.1 μL | 4671080BT |
| • E1-ClipTip Equalizer 8-ch 10-300 μL | 10-300 | 9.0-14.2 | 0.1 μL | 4672080BT |
| 300 μL Tips | | | Packaging | |
| • ClipTip 300 | 10-300 | | Racked, 96 tips/rack, 10 racks/pack | 94410510 |
| • ClipTip 300 | 10-300 | | Racked, sterile, 96 tips/rack, 10 racks/pack | 94410513 |
| • ClipTip 300 | 10-300 | | Reload tower, 96 tips/insert, 10 inserts/pack | 94410517 |
| • ClipTip 300 | 10-300 | | Reload tower, Sterile, 96 tips/insert, 10 inserts/pack | 94410518 |
| • ClipTip 300 | 10-300 | | Racked, filtered, sterile, 96 tips/rack, 10 racks/pack | 94420513 |
| • ClipTip 300 | 10-300 | | Reload insert, filtered, sterile, 96 tips/insert, 10 inserts/pack | 94420518 |
| • ClipTip 300 Ext | 10-300 | | Racked, 96 tips/rack, 8 racks/pack | 94410610 |
| • ClipTip 300 Ext | 10-300 | | Racked, sterile, 96 tips/rack, 8 racks/pack | 94410613 |
| • ClipTip 300 Ext | 10-300 | | Reload tower, 96 tips/insert, 10 inserts/pack | 94410617 |
| • ClipTip 300 Ext | 10-300 | | Racked, filtered, sterile, 96 tips/rack, 8 racks/pack | 94420613 |
| 850 μL Pipettes | | | Increment | |
| • E1-ClipTip 12-ch 30-850 μL | 30-850 | | 1.0 μL | 4671090 |
| 1000 μL Pipettes | | | Packaging | |
| • ClipTip 1000 | 30-1000 | | Racked, 96 tips/rack, 8 racks/pack | 94410710 |
| • ClipTip 1000 | 30-1000 | | Racked, sterile, 96 tips/rack, 8 racks/pack | 94410713 |
| • ClipTip 1000 | 30-1000 | | Reload tower, 96 tips/insert, 8 inserts/pack | 94410717 |
| • ClipTip 1000 | 30-1000 | | Reload tower, sterile, 96 tips/insert, 8 inserts/pack | 94410718 |
| • ClipTip 1000 | 30-1000 | | Racked, filtered, sterile, 96 tips/rack, 8 racks/pack | 94420713 |
| • ClipTip 1000 | 30-1000 | | Reload insert, filtered, sterile, 96 tips/insert, 8 inserts/pack | 94420718 |
| 1250 μL Pipettes | | | Increment | |
| • E1-ClipTip 15-1250 μL | 15-1250 | | 1.0 μL | 4670040BT |
| • E1-ClipTip 8-ch 15-1250 μL | 15-1250 | | 1.0 μL | 4671100BT |
| • E1-ClipTip Equalizer 6-ch 15-1250 μL | 15-1250 | 9.0-19.8 | 1.0 μL | 4672090BT |
| • E1-ClipTip Equalizer 8-ch 15-1250 μL | 15-1250 | 9.0-14.2 | 1.0 μL | 4672100BT |
| 1250 μL Tips | | | Packaging | |
| • ClipTip 1250 | 15-1250 | | Racked, 96 tips/rack, 8 racks/pack | 94410810 |
| • ClipTip 1250 | 15-1250 | | Racked, sterile, 96 tips/rack, 8 racks/pack | 94410813 |
| • ClipTip 1250 | 15-1250 | | Reload insert, 96 tips/insert, 8 inserts/pack | 94410817 |
| • ClipTip 1250 | 15-1250 | | Reload insert, sterile, 96 tips/insert, 8 inserts/pack | 94410818 |
| • ClipTip 1250 | 15-1250 | | Racked, filtered, sterile, 96 tips/rack, 8 racks/pack | 94420813 |
| • ClipTip 1250 | 15-1250 | | Reload insert, filtered, sterile, 96 tips/insert, 8 inserts/pack | 94420818 |

* ClipTip 384 and ClipTip 12.5 tips and compatible E1-ClipTip models utilize snap and clip tip fitting technology.

All ClipTip pipette tips are low retention with the exception of 384 and 12.5 models.

** Bluetooth models available in following countries: Andorra, Australia, Austria, Belgium, Bulgaria, Canada, China, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hong Kong, Hungary, Iceland, India, Ireland, Italy, Japan, Korea (Rep. of), Latvia, Liechtenstein, Lithuania, Luxembourg, Malta, Monaco, Netherlands, New Zealand, Norway, Poland, Portugal, Romania, San Marino, Singapore, Slovakia, Slovenia, Spain, Sweden, Switzerland, Taiwan, United Kingdom, United States and Vatican. For other countries catalog numbers without 'BT' extension should be used.


| | | |
|---|---|--|
| • 2 μL /10 μL / 12.5 μL / 20 μL | • 100 μL / 125 μL / 200 μL | • 850 μL / 1000 μL |
| • 30 μL / 50 μL | • 300 μL | • 1250 μL |

Revolutionize the way you pipette with a new web-based pipetting app

E1-ClipTip Electronic Pipettes and ClipTip Pipette Tips



E1-ClipTip multichannel pipette



E1-ClipTip Equalizer pipette with adjustable tip spacing
(See chart on the back for more details on adjustable tip spacing)



New CT300 Extended Tips
Enables reach to the bottom of test tubes



E1-ClipTip single charging stand



E1-ClipTip multi charging stand



Thermo Scientific Good Laboratory Pipetting

Learn how the right combination of technique, tools, ergonomics and service support Good Laboratory Pipetting at

thermofisher.com/pipette



Accessories for ClipTip Pipettes and Tips - ordering information









| Description | Color Key | Increment/Packaging | Cat. No. |
|---|-----------|-------------------------------------|----------|
| F-Series Stand, white, 6 positions | | | 9420400 |
| E1-ClipTip charging stand, 1 position | | Compatible with E1-ClipTip pipettes | 9420500 |
| E1-ClipTip multi charging stand, 3 position | | Compatible with E1-ClipTip pipettes | 9420510 |
| Empty rack, small | • • | 10 racks, no tips/inserts | 94410219 |
| Empty rack, medium | • • | 10 racks, no tips/inserts | 94410519 |
| 300 Ext empty rack | • | 8 racks, no tips/inserts | 94410619 |
| Empty rack, large | • • | 8 racks, no tips/inserts | 94410819 |

• 2 μ L / 10 μ L / 12.5 μ L / 20 μ L
 • 100 μ L / 125 μ L / 200 μ L
 • 850 μ L / 1000 μ L
• 30 μ L / 50 μ L
 • 300 μ L
 • 1250 μ L

Save 90% more time with E1-ClipTip Equalizer pipettes!

Save over 90% more time, and 90% less repetitions. Transfer samples between any labware format, and save time with E1-ClipTip electronic adjustable tip spacing models.

Equalizer model selector guide

| 96-Format Equalizer Pipettes | | | | | |
|-------------------------------|--------------|---|---|--|--------------|
| Channels | Max. spacing | Transfer to/from | ↔ | Transfer from/to | Min. spacing |
| 6-channel | 19.8 mm | <ul style="list-style-type: none"> Test tube racks 24-well microplates  | ↔ | <ul style="list-style-type: none"> 96-well microplates Deepwell blocks  | 9.0 mm |
| 8-channel | 14.2 mm | <ul style="list-style-type: none"> Microcentrifuge tube racks 48-well microplates  | ↔ | <ul style="list-style-type: none"> 96-well microplates Deepwell blocks  | 9.0 mm |
| 384-Format Equalizer Pipettes | | | | | |
| Channels | Max. spacing | Transfer to/from | ↔ | Transfer from/to | Min. spacing |
| 8-channel | 14.2 mm | <ul style="list-style-type: none"> 96-well microplates 48-well microplates Deepwell blocks Microcentrifuge tube racks  | ↔ | <ul style="list-style-type: none"> 384-well microplates 384 PCR plates Agarose gels  | 4.5 mm |
| 12-channel | 9.0 mm | <ul style="list-style-type: none"> 96-well microplates Deepwell blocks  | ↔ | <ul style="list-style-type: none"> 384-well microplates 384 PCR plates Agarose gels  | 4.5 mm |



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.