



Ultrawave

Precision ultrasonic cleaning equipment

QS13x

13 Litre Ultrasonic Bath

Ultrawave's **QS13x** offers outstanding performance for use in scientific and laboratory applications, medical, chiropody, veterinary and many other professional environments.

Incorporating Ultrawave's *Frequency Leap* Technology, homogenous cavitation activity is generated throughout the ultrasonic bath reducing dead spots and standing waves. This provides a more consistent scrubbing action for excellent results.



The **QS13x** is a 13 Litre heated, digital ultrasonic bath.

The menu driven control panel allows simple adjustment of cycle time (0 to 99 minutes), temperature (ambient to 70°C), ultrasonic power and frequency mode for precise and accurate processing.

Combining features such as degas mode and fluid level sensor, the **QS13x** offers repeatable cleaning performance with consistent cleaning results every time.



QS13x

13 Litre Ultrasonic Bath

Construction Case	304 Stainless Steel
Construction Tank	304 Stainless Steel
Front Panel	ABS Plastic
Control	Digital LCD
Overall Dimensions	375 x 345 x 265mm (WxDxH)
Tank Dimensions	330 x 300 x 150mm (WxDxH)
Basket Dimensions	285 x 260 x 110mm (WxDxH)
Maximum Capacity	13.0Litres
Working Capacity	10.5Litres
Ultrasonic Power	200W
Ultrasonic Power (per Litre)	19W/Litre
Operating Frequency	32 - 38KHz <i>Frequency Leap</i>
Maximum Heating Capacity	Ambient - 70°C
Heater Power	300W
Timer	0 - 99 minutes
Mains Connection	230V(50Hz)
Drain Outlet	½" BSP VALVE
Lid	ABS Plastic
Basket	Stainless Steel
Gross Weight (inc. lid, basket & packaging)	10.9kg
Net Weight	8.25kg
Warranty	3 Years





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.