

# Qi400

## 30Litre Ultrasonic Bath

Ultrawave's **Qi-Series** are robust and reliable, delivering outstanding cleaning performance in a wide range of industries. From precision, pharmaceutical and medical device manufacturing to heavy duty cleaning in automotive and engineering processing. **Qi-Series** Ultrasonic Cleaners are Class I Medical Devices, compliant to MDR (EU) 2017/745.

The 2mm thick stainless steel tank provides longevity and durability whilst excellent ultrasonic activity is generated via Ultrawave's *Frequency Leap* technology.

*Frequency Leap* technology reduces dead spots and promotes homogenous cavitation activity within the tank. This provides a more consistent scrubbing action for highly effective cleaning results.



The **Qi-400** is a 30Litre ultrasonic cleaner with a menu driven, LCD control panel. This allows simple adjustment of cycle time (0 to 99mins or infinity mode), temperature (ambient to 80°C), ultrasonic power and frequency mode for precise and accurate processing.

Incorporating features such as degas mode and a fluid level sensor, the

**Qi-Series** is reliable and durable, giving you consistent cleaning results every time.

Specifications are subject to change for continual improvements

# Qi400

## 30Litre Ultrasonic Bath

Construction Case	304 Stainless Steel
Construction Tank	316 Stainless Steel
Control	Digital LCD
Overall Dimensions	530 x 380 x 395mm (WxDxH)
Tank Dimensions	485 x 275 x 220mm (WxDxH)
Internal Basket Dimensions	455 x 250 x 140mm (WxDxH)
Maximum Capacity	30Litres
Working Capacity	26Litres
Ultrasonic Power	400W
Ultrasonic Power (per litre)	15W/Litre
Operating Frequency	32 - 38kHz <i>Frequency Leap</i>
Maximum Heating Capacity	Ambient - 80°C
Heater Power	1350W
Timer	0 - 99 minutes or Continuous
Mains Connection	230VAC
Drain Outlet	Lockable 1/2" BSP valve (located at rear of machine)
Lid	304 Stainless Steel
Basket	304 Stainless Steel
Gross Weight (inc. lid, basket & packaging)	24.5kg
Net Weight	22kg
Warranty	1 Year as standard - extend to 3 Years if registered online with us (within 30days of purchase)



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.