

U50

0.5Litre Ultrasonic Bath

Ultrawave's **U-Series** is compact and user friendly, offering effective cleaning and processing of small items in both professional and domestic applications.

These include healthcare, podiatry, tattooists, jewellers and many more.

Available in heated and unheated versions, the **U-Series** allows simple control of time (0 to 15 minutes) and temperature (ambient to 70°C) where heated options are selected.



The **U50** is a 0.5Litre ultrasonic bath with a timer control of up to 15 minutes.

Each **U-Series** is supplied complete with a stainless steel basket and ABS plastic lid, useable as a drip tray to collect excess water when the basket is removed from the tank.

The **U50** provides the cost benefits of excellent entry level ultrasonics with Ultrawave's long renowned quality and reliability.



U50

0.5Litre Ultrasonic Bath

Construction Case	304 Stainless Steel
Construction Tank	304 Stainless Steel
Control	Analogue (Time)
Overall Dimensions	160 x 95 x 165mm (WxDxH)
Tank Dimensions	150 x 85 x 65mm (WxDxH)
Basket Dimensions	120 x 60 x 25mm (WxDxH) (50mm height with handles)
Maximum Capacity	0.6Litres
Working Capacity	0.5Litres
Ultrasonic Power	35W
Ultrasonic Power (per litre)	70W/litre
Operating Frequency	44KHz
Maximum Heating Capacity	N/A
Heater Power	N/A
Timer	0 - 15 minutes
Mains Connection	220 – 240V (50 - 60Hz)
Drain Outlet	N/A
Lid	ABS Plastic
Basket	Stainless Steel
Gross Weight (inc. lid, basket & packaging)	1.8KGs
Net Weight	1.4KGs
Warranty	3 Years



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.