

# SIDE ACCESS RACKS

Side Opening Horizontal units also known as side access racks, fit onto shelves of upright freezers to hold cryoboxes or microtiter plates. Side access racks feature a handle on one end to remove it from the freezer, these racks are best suited for longer term storage. Side access racks are manufactured from aluminium or stainless steel. Side access racks can be any size to meet the requirements of the user. They can also feature a locking rod to hold items in place, and have a handle fitted to the top to make the unit easier to transport.

Dimensions (mm)	Box Size / Boxes Per Rack	Part Number
295H x 145W x 290D	4" (100mm) / 4	<b>SOR455</b>
295H x 145W x 290D	3" (75mm) / 6	<b>SOR453</b>
295H x 145W x 290D	2" (50mm) / 10	<b>SOR110</b>
420H x 140W x 560D	4" (100mm) / 12	<b>SOR304</b>
420H x 140W x 560D	3" (75mm) / 20	<b>SOR368</b>
420H x 140W x 560D	2" (50mm) / 28	<b>SOR090</b>
235H x 141W x 560D	2" (50mm) / 16	<b>SOR196</b>
235H x 141W x 560D	3" (75mm) / 12	<b>SOR197</b>
290H x 140W x 560D	2" (50mm) / 20	<b>SOR032</b>
340H x 140W x 560D	2" (50mm) / 24	<b>SOR033</b>
230H x 140W x 560D	2" (50mm) / 16	<b>SOR196</b>

*Also available in stainless steel.*


Wesbart also offers a very popular adjustable side access rack which offers the user flexibility to store different size boxes. Using the unique design of this rack the user can set the removable shelves to store either 16x 2" boxes, 12x 3" boxes or 8x 4" boxes. *Part Number SOR002*

## BESPOKE RACKS


Should you require a bespoke front or side access rack please visit our website ([www.wesbart.co.uk](http://www.wesbart.co.uk)) to fill in a quotation request. The information that will be required is the freezer make and model, maximum internal size of the shelf and the size of the item you require racking.



Side access rack to store 2" box **SOR196**



Side access rack to store microtiter plates **SOR004**



Side access rack to store various height boxes **SOR002B**



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.